

### What You Need to Do

1. Pure cultures should be supplied inoculated onto solid media and incubated to produce sufficient growth.
2. Please seal slopes or plates with parafilm and then wrap in bubble wrap and place in a plastic tub inside a cardboard box. Send by courier to the following address:

nadicom GmbH  
Hertzstr. 16, Gebäude 6.40  
D-76187 Karlsruhe  
Germany

Tel: (+49) 64 21 - 131 75 (Dr. Nüßlein)  
Tel. (+49) 721-6084 - 4627 (laboratory)  
Email: [info@nadicom.com](mailto:info@nadicom.com)

3. Please note that we can store cultures in glycerol stocks at -80°C for up to 3 years if requested.
4. Please complete the sequencing request form and contact name with each identification.



